




Quick Reference Guide


 This Quick Reference Guide does not replace the comprehensive User's Manual for your Accu-Chek Go meter. Please be sure to read the safety notices contained in the Instructions for use.

1.



2.






- Top and sides**
- ① button: turn on or off, start a set-up to change settings, save entries
 - ② Liquid crystal display (LCD)
 - ③ and buttons: change settings (decrement and increment value, respectively), retrieve results, request averages
 - ④ Test strip guide – this is where you insert the test strip
 - ⑤ Measurement optics cover
 - ⑥ Test strip ejector
 - ⑦ Code chip
 - ⑧ Infrared port for downloading results to a PC, for instance




- Bottom**
- ⑦ Code chip
 - ⑨ Code chip slot
 - ⑩ Type plate
 - ⑪ Battery compartment
 - ⑫ Battery (type CR 2032)



- 1 Set-up
- 2 When an average is displayed: number of days covered
- 3 LH or LL: target range setting (set-up)
- 4 Alarm clock function turned on. Flashing: optical signal during alarm
- 5 Test or averaging in progress
- 6 Unit of measurement for BG result
- 7 Flag indicating special result
- 8 Flashing: cue to allow test strip to draw up blood
- 9 Flashing: cue to insert test strip
- 10 Test strips are past their expiry date
- 11 Code number display. Flashing: code chip is missing
- 12 Target range setting (set-up). Flashing: result is outside the target range
- 13 Flag indicating quality control
- 14 Flashing: temperature is/was outside the permitted range
- 15 Beep tone or Acoustic Mode is turned on
- 16 Flashing: battery is running low
- 17 When an average is displayed: number of results included
- 18 Memory, averages or results download
- 19 Follows time in 12-hour format

- 20 Display messages**
- Code chip is missing. In memory: no results in memory, averaging not possible or result erased
 - AL Alarm setting (set-up)
 - CLr Flashing: latest result being erased
 - End End of set-up, self-test following test or download complete
 - F In Quick set-up: group selected
 - H_i Result >600 mg/dL (33.3 mmol/L)
 - L_o Result <10 mg/dL (0.6 mmol/L)
 - PC Results download

Quick set-up: Changing settings



- Error messages**
- E-1 Test strip inserted incorrectly, test strip bent/moved, strong electromagnetic field
 - E-2 Code chip from another system, meter cannot read code chip, code chip defective, code chip removed during test
 - E-4 Measurement optics dirty, test strip drew up blood too soon, test strip bent/moved, used strip inserted, test strip past its expiry date, new test strip inserted too soon after last test (while End was showing)
 - E-6 Test strip drew up too little blood, test strip drew up blood too soon, test strip bent/moved, battery empty, meter error, strong electromagnetic field
 - EEE Meter error
 - ttt Temperature too low or too high to test.